Trephine Ejection Kit

DE

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Trephine Ejection Kit

EN

Trephine System for Extraction of Cylindrical Bone Grafts
Developed by Prof. Dr. Fouad Khoury

Art.-No.: STE00

The Trephine Ejection Kit was developed together with Prof. Dr. Fouad Khoury. It contains several unique ideas which turn this kit into something special. Two-part, internally cooled trepanes allow for a safe and minimal invasive extraction of local bone cylinders. The inner cooling helps to prevent a damage of the bone tissue due to overheating, while the removable trephine working parts allow for an especially easy sampling of the cut out bone cylinders. Either with the aid of the ejection needle or the ejection instruments, which perfectly match with the diameter of the trephine working parts, the bone cylinders can be pushed out of the trepanes. Furthermore, the included pre-drilling trepanes allow for a safe application of the trepanes. Due to the two-part trepanes, the variety of different diameters and the optimal matching of the instruments, the Trephine Ejection Kit offers highest efficiency, flexibility and precision to the user for the extraction of bone cylinders.

Indications

The Bone Management® system Trephine Ejection Kit is indicated for the intraoral extraction of bone cylinders for augmentation of bone deficits in the upper and lower jaw.

Instruction

The recommended speeds for the application and the maximum speeds are summarized in the instrument overview under „Content“. For proper, safe use, these must be observed.

Before Application

The instruments are supplied non-sterile, so they need to be cleaned, disinfected and sterilized prior to use, this applies in particular to the first-time use. For cleaning and disinfection the products are disassembled as far as possible. For that the shank of the trephine is attached to the hand wheel, so the trephine working part can be released and removed by a left-hand turn. The rubber ring remains on the shank during the first processing. After clearing and disinfection have been completed, the trephine working part is attached again to the shank and fixed by a right-hand turn. The secure fit of the working part has to be checked. Subsequently, the instruments can be sterilized according to the instructions for processing. Since the trepanes have an internal cooling, these are instruments with increased processing requirements. According to the RIKI guideline, mechanical cleaning and disinfection in disinfection devices and sterilization with moist heat is recommended. For proper processing of the articles, please note the instructions for processing (cleaning, disinfection and sterilization) of medical devices produced by Hager & Meisinger GmbH.

Initial Drilling

According to the planned diameter of the bone cylinder the appropriate initial drill is selected. With this drill the bone at the donor site is punch-marked in order to avoid slippage of the trephine.

Extraction of the bone cylinder from the donor site

According to the planned diameter of the bone cylinder the appropriate trephine is selected. Before use, the fixed position of the trephine working part on the shank has to be checked.

Caution: Depending on the construction, the trepanes may not be used counter-clockwise. At the punch-marked donor site the bone cylinder of the planned length can now be extracted with the aid of the trephine. For this purpose, the trephine is placed on the cortical bone and intermittently inserted into the bone. It is recommended to support the working hand to control the trephine optimally. The speed should not be too low to avoid tilting of the trephine in the depth of the bone. Please note the recommended and maximum speeds. The depth marks on the trephine working parts serve for length orientation.

Extraction of the bone cylinder from the trephine
To extract the bone cylinder from the trephine working part, the trephine is removed from the handpiece. To achieve a better handling, the ejection needle is attached to the hand wheel. Then, it can be fed through the lumen of the trepan's internal cooling from the shank end in the direction of the trephine working part, so the bone cylinder is pushed out.

Caution: In the case of soft bone, the application of the ejection needle is contraindicated due to its small diameter. In this case, proceed according to step "Alternative Extraction of the bone cylinder from the trephine in case of soft bone".

Alternative Extraction of the bone cylinder from the trephine in case of soft bone

In case of soft bone, the extraction of the bone cylinder is performed with the aid of the ejection instruments. For that, the shank of the trephine is attached to the hand wheel, so that the trephine working part can be released and removed by a left-hand turn. Then, the trephine working part is inserted into the ejection sleeve. With the aid of the appropriate ejection instrument the bone cylinder is pushed out. The 2.1 mm and 2.5 mm diameter trepanes go with the 2.1 mm diameter ejection instrument. The 2.9 mm and 3.1 mm diameter trepanes go with the 2.9 mm diameter ejection instrument.

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After use

The instruments need to be cleaned, disinfected and sterilized after every application. For cleaning and disinfection the products are disassembled as far as possible...

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Contraindications

Basically, general medical as well as local, absolute and relative contraindications for dental surgical procedures must be considered.

Absolute contraindications

- Dentoalveolar growth which has not come to an end (exception: cases in which no dentoalveolar growth can be expected, e.g. ectodermal dysplasia)
Active infections as well as pathological processes
Insufficient bone supply (quality / quantity) of the donor site as well as insufficient bone quality of the donor region

Relative contraindications

- Diseases that affect bone metabolism
Drug and alcohol abuse
Lack of cooperation of the patient
Poor circulation
Hard work or active sports
Mental condition, which can lead to disregard of the medical order
Highly atrophic jaw

IMPORTANT

Attention must be paid to protecting anatomical structures (safety clearance at least 2 mm) as well as the gradient of the adjacent teeth / tooth roots (risk of damage, infection / dehiscences).

CAUTION

For all cutting and drilling procedures (transport removal / other work), the following applies: In order to reduce the risk of overheating the bone and necrosis formation along with it, the respective instrument has to be used intermittently using light pressure and continuous cooling using a sterile physiological saline solution.

General instructions

Please follow general application and safety instructions for MEISINGER products in the medical area and also the advice for processing (cleaning, disinfection and sterilization) of medical devices from Hager & Meisinger GmbH.

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Logo Hager & Meisinger GmbH, CE 0044, contact information for Germany.