|     |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
|-----|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| en  | <b>Instructions for use/Technical description</b><br>Multi-part instruments                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| USA | <b>Note for U.S. users</b><br>This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at <a href="http://www.aesculapusaifus.com">www.aesculapusaifus.com</a> . If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost. |
| de  | <b>Gebrauchsanweisung/Technische Beschreibung</b><br>Mehrteilige Instrumente                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| fr  | <b>Mode d'emploi/Description technique</b><br>Instruments composés de plusieurs pièces                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| es  | <b>Instrucciones de manejo/Descripción técnica</b><br>Instrumentos compuestos                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| it  | <b>Istruzioni per l'uso/Descrizione tecnica</b><br>Strumenti multiparte                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| pt  | <b>Instruções de utilização/Descrição técnica</b><br>Instrumentos de várias peças                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| nl  | <b>Gebruiksaanwijzing/Technische beschrijving</b><br>Meerdelige instrumenten                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| da  | <b>Brugsanvisning/Teknisk beskrivelse</b><br>Flerdelte instrumenter                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| sv  | <b>Bruksanvisning/Teknisk beskrivning</b><br>Flerdelade instrument                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| fi  | <b>Käyttöohje/Tekninen kuvaus</b><br>Moniosaisset instrumentit                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| lv  | <b>Lietošanas instrukcijas/tehniskais apraksts</b><br>Daudzdaļīgi instrumenti                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| lt  | <b>Naudojimo instrukcija/techninis aprašas</b><br>Sudėtiniai instrumentai                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| ru  | <b>Инструкция по применению/Техническое описание</b><br>Многокомпонентные инструменты                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| cs  | <b>Návod k použití/Technický popis</b><br>Více dílné nástroje                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| pl  | <b>Instrukcja użytkowania/Opis techniczny</b><br>Instrumenty wieloczęściowe                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| sk  | <b>Návod na použitie/Technický opis</b><br>Viacdielne nástroje                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| hu  | <b>Használati útmutató/Műszaki leírás</b><br>Többrészes műszerek                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| sl  | <b>Navodila za uporabo/Tehnični opis</b><br>večdelnimi instrumenti                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| hr  | <b>Upute za uporabu/Tehnički opis</b><br>Višedijelni instrumenti                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| ro  | <b>Manual de utilizare/Descriere tehnică</b><br>Instrumente multipartite                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| bg  | <b>Упътване за употреба/Техническо описание</b><br>Многокомпонентни инструменти                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| tr  | <b>Kullanım Kılavuzu/Teknik açıklama</b><br>Çok parçalı aletler                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| el  | <b>Οδηγίες χρήσης/Τεχνική περιγραφή</b><br>Εργαλεία πολλαπλών μερών                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany  
Phone +49 (0) 7461 95-0 | Fax +49 (0) 7461 95-26 00 | [www.bbraun.com](http://www.bbraun.com)

AESCULAP® – a B. Braun brand

TA013788    2021-04    Change No. AE0060622

1. About this document

*Note*  
General risk factors associated with surgical procedures are not described in these instructions for use.

1.1 Scope

Multi-part instruments, all surfaces of which are directly accessible and visible after intuitive disassembly (e.g. screwing, inserting, pulling). The instruments do not contain aluminium, nickel-silver or tin components and are not nickel-plated or chrome-plated.

*Note*  
The applicable CE mark for the product can be seen on the label or packaging of the product.

- For article specific instructions for use and material compatibility and lifetime information, see B. Braun eIFU at eifu.bbraun.com

1.2 Safety messages

Safety messages make clear the dangers to patient, user and/or product that could arise during the use of the product. Safety messages are labeled as follows:

**⚠ WARNING**  
Indicates a possible threat of danger. If not avoided, minor or moderate injury may result.

**⚠ CAUTION**  
Indicates a possible threat of material damage. If not avoided, the product may be damaged.

2. Clinical use

2.1 Areas of use and limitations of use

2.1.1 Intended use

The instruments are intended for universal surgical use in all disciplines.

2.1.2 Indications

*Note*  
The manufacturer is not responsible for any use of the product against the specified indications and/or the described applications.  
  
For indications, see Intended use.

2.1.3 Absolute contraindications

No known absolute contraindications.

2.1.4 Relative contraindications

The following conditions, individual or combined, can lead to delayed healing or compromise the success of the operation:  
■ Medical or surgical conditions (e.g. comorbidities) which could hinder the success of the operation.  
In the presence of relative contraindications, the user decides individually regarding the use of the product.

2.2 Safety information

2.2.1 Clinical user

**General safety information**  
To prevent damage caused by improper setup or operation, and to not compromise the manufacturer warranty and liability:  
► Use the product only according to these instructions for use.  
► Follow the safety and maintenance instructions.  
► Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge and experience.  
► Store any new or unused products in a dry, clean, and safe place.  
► Prior to use, check that the product is in good working order.  
► Keep the instructions for use accessible for the user.

*Note*  
The user is obligated to report all severe events in connection with the product to the manufacturer and the responsible authorities of the state in which the user is located.

**Notes on surgical procedures**  
It is the user's responsibility to ensure that the surgical procedure is performed correctly.  
Appropriate clinical training as well as a theoretical and practical proficiency of all the required operating techniques, including the use of this product, are prerequisites for the successful use of this product.  
The user is required to obtain information from the manufacturer if there is an unclear preoperative situation regarding the use of the product.

2.2.2 Sterility

The product is delivered in an unsterile condition.  
► Clean the new product after removing its transport packaging and prior to its initial sterilization.

2.3 Application

**⚠ WARNING**  
**Risk of injury and/or malfunction!**  
► Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.  
► Always carry out a function test prior to each use of the product.

3. Validated reprocessing procedure

3.1 General safety information

*Note*  
Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.  
  
*Note*  
For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

*Note*  
Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

*Note*  
Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.

*Note*  
If there is no final sterilization, then a virucidal disinfectant must be used.

*Note*  
For up-to-date information about reprocessing and material compatibility, see B. Braun eIFU at eifu.bbraun.com  
The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

3.2 General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.  
Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel.  
Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying.  
Additional drying, if necessary.  
Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:  
■ Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.  
■ Material damage such as corrosion, cracks, fracturing, premature aging or swelling.  
► Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.  
► Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.a-k-i.org, link to "AKI-Brochures", "Red brochure".

3.3 Reusable products

Influences of the reprocessing which lead to damage to the product are not known.  
A careful visual and functional inspection before the next use is the best opportunity to recognize a product that is no longer functional, see Inspection.

3.4 Preparations at the place of use

- If applicable, rinse non-visible surfaces preferably with deionized water, with a disposable syringe for example.
- Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

3.5 Preparing for cleaning

- Disassemble the product prior to cleaning.

3.6 Cleaning/Disinfection

**3.6.1 Product-specific safety information on the reprocessing method**  
Damage to or destruction of the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!  
► Use cleaning agents and disinfectants according to the manufacturer's instructions.  
► Observe specifications regarding concentration, temperature and exposure time.  
► Do not exceed the maximum allowable disinfection temperature of 95 °C.

3.6.2 Validated cleaning and disinfection procedure

| Validated procedure                                                                                 | Specific requirements                                                                                                                                                                          | Reference                                                                                                                                                                                          |
|-----------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Manual cleaning with immersion disinfection                                                         | ■ Suitable cleaning brush<br>■ Disposable syringe 20 ml<br>■ Drying phase: Use a lint-free cloth or medical compressed air                                                                     | Chapter Manual cleaning/disinfection and subsection:<br>■ Chapter Manual cleaning with immersion disinfection                                                                                      |
| Manual pre-cleaning with brush and subsequent mechanical alkaline cleaning and thermal disinfection | ■ Suitable cleaning brush<br>■ Disposable syringe 20 ml<br>■ Mount jaws protection on the product.<br>■ Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots). | Chapter Mechanical cleaning/disinfection with manual pre-cleaning and subsection:<br>■ Chapter Manual pre-cleaning with a brush<br>■ Chapter Mechanical alkaline cleaning and thermal disinfecting |

3.7 Manual cleaning/disinfection

- Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.
- After manual cleaning/disinfection, check visible surfaces visually for residues.
- Repeat the cleaning/disinfection process if necessary.

3.7.1 Manual cleaning with immersion disinfection

| Phase | Step                  | T [°C/°F] | t [min] | Conc. [%] | Water quality | Chemical                                                       |
|-------|-----------------------|-----------|---------|-----------|---------------|----------------------------------------------------------------|
| I     | Disinfecting cleaning | RT (cold) | >15     | 2         | D-W           | Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9* |
| II    | Intermediate rinse    | RT (cold) | 1       | -         | D-W           | -                                                              |
| III   | Disinfection          | RT (cold) | 5       | 2         | D-W           | Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9* |
| IV    | Final rinse           | RT (cold) | 1       | -         | FD-W          | -                                                              |
| V     | Drying                | RT        | -       | -         | -             | -                                                              |

D-W: Drinking water  
FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)  
RT: Room temperature  
\*Recommended: BBraun Stabimed fresh  
  
► Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

**Phase I**  
► Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.  
► Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.  
► If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.  
► Mobilize non-rigid components, such as set screws, links, etc. during cleaning.  
► Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

- Phase II**
- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
  - ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
  - ▶ Drain any remaining water fully.
- Phase III**
- ▶ Fully immerse the product in the disinfectant solution.
  - ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
  - ▶ Rinse lumens at least 5 times at the beginning of the exposure time using an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.
- Phase IV**
- ▶ Rinse/flush the product thoroughly (all accessible surfaces).
  - ▶ Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
  - ▶ Rinse lumens with an appropriate disposable syringe at least five times.
  - ▶ Drain any remaining water fully.
- Phase V**
- ▶ Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

3.8 Mechanical cleaning/disinfection with manual pre-cleaning

*Note*  
The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

*Note*  
The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

3.8.1 Manual pre-cleaning with a brush

| Phase | Step                  | T<br>[°C/°F] | t<br>[min] | Conc.<br>[%] | Water<br>quality | Chemical                                                       |
|-------|-----------------------|--------------|------------|--------------|------------------|----------------------------------------------------------------|
| I     | Disinfectant cleaning | RT<br>(cold) | >15        | 2            | D-W              | Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9* |
| II    | Rinsing               | RT<br>(cold) | 1          | -            | D-W              | -                                                              |

D-W: Drinking water  
RT: Room temperature  
\*Recommended: BBraun Stabimed fresh

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

- Phase I**
- ▶ Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
  - ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
  - ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
  - ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
  - ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

- Phase II**
- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
  - ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

3.8.2 Mechanical alkaline cleaning and thermal disinfecting  
Machine type: single-chamber cleaning/disinfection device without ultrasound

| Phase | Step                 | T<br>[°C/°F] | t<br>[min] | Water<br>quality | Chemical                                                                                                                                                                                                                               |
|-------|----------------------|--------------|------------|------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| I     | Prerinse             | <25/77       | 3          | D-W              | -                                                                                                                                                                                                                                      |
| II    | Cleaning             | 55/131       | 10         | FD-W             | <div><div>■</div>Concentrate, alkaline:<ul style="list-style-type: none"><li>- pH ~ 13</li><li>- &lt;5 % anionic surfactant</li></ul><div>■</div>working solution 0.5%<ul style="list-style-type: none"><li>- pH = 11*</li></ul></div> |
| III   | Intermediate rinse   | >10/50       | 1          | FD-W             | -                                                                                                                                                                                                                                      |
| IV    | Thermal disinfecting | 90/194       | 5          | FD-W             | -                                                                                                                                                                                                                                      |
| V     | Drying               | -            | -          | -                | In accordance with the program for the cleaning and disinfecting machine                                                                                                                                                               |

D-W: Drinking water  
FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

\*Recommended: BBraun Helimatic Cleaner alkaline

- ▶ Check visible surfaces for residues after mechanical cleaning/disinfecting.

3.9 Inspection

- ▶ Allow the product to cool down to room temperature.
- ▶ Dry the product if it is wet or damp.

3.9.1 Visual inspection

- ▶ Ensure that all soiling has been removed. In particular, pay attention to mating surfaces, hinges, shafts, recessed areas, drill grooves and the sides of the teeth on rasps.
- ▶ If the product is dirty: repeat the cleaning and disinfection process.
- ▶ Check the product for damage, e.g. insulation or corroded, loose, bent, broken, cracked, worn or severely scratched and fractured components.
- ▶ Check the product for missing or faded labels.
- ▶ Check the products with long, slim shapes (in particular rotating instruments) for deformities.
- ▶ Check the product for damage to the spiral element.
- ▶ Check the cutting edges for continuity, sharpness, nicks and other damage.
- ▶ Check the surfaces for rough spots.
- ▶ Check the product for burrs that could damage tissue or surgical gloves.
- ▶ Check the product for loose or missing parts.
- ▶ Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical service.

3.9.2 Functional test

- ⚠ CAUTION**  
Damage (metal cold welding/friction corrosion) to the product caused by insufficient lubrication!
- ▶ Oil moving parts (e.g. joints, spool parts and threaded rods) prior to the functional test using oil suitable for the applied sterilization process (e.g. STERILIT® I oil spray JG600 for steam sterilization or STERILIT® I drip-feed lubricator JG598).
  - ▶ Check that the product functions correctly.
  - ▶ Check that all moving parts are working property (e.g. hinges, locks/latches, sliding parts etc.).
  - ▶ Check for compatibility with associated products.
  - ▶ Immediately put aside inoperative products and send them to Aesculap Technical Service, see Technical service.

3.10 Packaging

- ▶ Appropriately protect products with fine working tips.
- ▶ Place the product in its holder or on a suitable tray. Ensure that sharp edges are covered.
- ▶ Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- ▶ Ensure that the packaging provides sufficient protection against contamination of the product during storage.

3.11 Steam sterilization

- Note*  
The product can be sterilized either in disassembled or in assembled condition.
- ▶ Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g., by opening any valves and faucets).
  - ▶ Validated sterilization process
    - Steam sterilization in fractionated vacuum process
    - Steam sterilizer in accordance with DIN EN 285 and validated in accordance with DIN EN ISO 17665
    - Sterilization in fractionated vacuum process at 134 °C, holding time 5 min
  - ▶ If several devices are sterilized at the same time in the same steam sterilizer: Ensure that the maximum permitted load according to the manufacturers' specifications is not exceeded.

3.12 Storage

- ▶ Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

4. Technical service

- ⚠ CAUTION**  
Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.
- ▶ Do not modify the product.
  - ▶ For service and repairs, please contact your national B. Braun/Aesculap agency.

**Service addresses**  
Aesculap Technischer Service  
Am Aesculap-Platz  
78532 Tuttlingen / Germany  
Phone: +49 7461 95-1601  
Fax: +49 7461 16-2887  
E-Mail: ats@aesculap.de  
Other service addresses can be obtained from the address indicated above.

5. Disposal

- ⚠ WARNING**  
Risk of infection due to contaminated products!
- ▶ Adhere to national regulations when disposing of or recycling the product, its components and its packaging.

- ⚠ WARNING**  
Risk of injury due to sharp-edged and/or pointed products!
- ▶ When disposing of or recycling the product, ensure that the packaging prevents injury by the product.

*Note*  
The user institution is obliged to reprocess the product before its disposal, see Validated reprocessing procedure.

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|     |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|-----|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| en  | <b>Instructions for use/Technical description</b><br>Tubular stainless steel instruments                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| USA | <b>Note for U.S. users</b><br>This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at <a href="http://www.aesculapusa.com">www.aesculapusa.com</a> . If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost. |
| de  | <b>Gebruiksaanweisung/Technische Beschreibung</b><br>Rohrförmige Instrumente aus Edelstahl                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| fr  | <b>Mode d'emploi/Description technique</b><br>Instruments de forme tubulaire en acier inoxydable                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| es  | <b>Instrucciones de manejo/Descripción técnica</b><br>Instrumental de acero inoxidable en forma de tubo                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| it  | <b>Istruzioni per l'uso/Descrizione tecnica</b><br>Strumenti tubolari in acciaio inossidabile                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| pt  | <b>Instruções de utilização/Descrição técnica</b><br>Instrumentos tubulares de aço inoxidável                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| nl  | <b>Gebruiksaanwijzing/Technische beschrijving</b><br>Buisvormige instrumenten van roestvrij staal                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| da  | <b>Brugsanvisning/Teknisk beskrivelse</b><br>Rørformede instrumenter af rustfrit stål                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| sv  | <b>Bruksanvisning/Teknisk beskrivning</b><br>Rörformade instrument i rostfritt stål                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| fi  | <b>Käyttöohje/Tekninen kuvaus</b><br>Putkimaiset instrumentit ruostumatonta terästä                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| lv  | <b>Lietošanas instrukcijas/tehniskais apraksts</b><br>Cauruļu formas instrumenti no nerūsējošā tērauda                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| lt  | <b>Naudojimo instrukcija/techninis aprašas</b><br>Vamzdeliniai nerūdijančio plieno instrumentai                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| ru  | <b>Инструкция по применению/Техническое описание</b><br>Трубчатые инструменты из высококачественной стали                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| cs  | <b>Návod k použití/Technický popis</b><br>Trubkové nástroje z ušlechtilé oceli                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| pl  | <b>Instrukcja użytkowania/Opis techniczny</b><br>Instrumenty rurkowe wykonane ze stali szlachetnej                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| sk  | <b>Návod na použitie/Technický opis</b><br>Trubicové nástroje z ušľachtilej ocele                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| hu  | <b>Használati útmutató/Műszaki leírás</b><br>Cső alakú eszközök rozsdamentes acélból                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| sl  | <b>Navodila za uporabo/Tehnični opis</b><br>Instrumenti v obliki cevi iz nerjavečega jekla                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| hr  | <b>Upute za uporabu/Tehnički opis</b><br>Cjevasti instrumenti od oplemenjenog čelika                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| ro  | <b>Manual de utilizare/Descriere tehnică</b><br>Instrumente tubulare din oțel inoxidabil                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| bg  | <b>Упътване за употреба/Техническо описание</b><br>Инструменти с форма на тръба, изработени от неръждаема стомана                                                                                                                                                                                                                                                                                                                                                                                                                           |
| tr  | <b>Kullanım Kılavuzu/Teknik açıklama</b><br>Paslanmaz çelikten yapılmış boru şeklindeki ekipmanlar                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| el  | <b>Οδηγίες χρήσης/Τεχνική περιγραφή</b><br>Σωληνωειδή εργαλεία από ανοξείδωτο χάλυβα                                                                                                                                                                                                                                                                                                                                                                                                                                                        |

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1.1 Scope

Tubular instruments such as aspirators and suction and irrigation cannulae that contain no aluminum, nickel silver or tin components and which are not nickel-plated or chrome-plated.

*Note*  
The applicable CE mark for the product can be found on the label or packaging of the product.

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2.1 Areas of use and limitations of use

2.1.1 Intended use

The surgical instruments are intended for universal use in various surgical disciplines.

2.1.2 Indications

*Note*  
The manufacturer is not responsible for any use of the product against the specified indications and/or the described applications.

For indications, see Intended use.

2.1.3 Contraindications

No known contraindications.

2.2 Safety information

2.2.1 Clinical user

General safety information

To prevent damage caused by improper setup or operation, and to not compromise the manufacturer warranty and liability:

- Use the product only according to these instructions for use.
- Follow the safety and maintenance instructions.
- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge and experience.
- Store any new or unused products in a dry, clean, and safe place.
- Prior to use, check that the product is in good working order.
- Keep the instructions for use accessible for the user.

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Notes on surgical procedures

It is the user's responsibility to ensure that the surgical procedure is performed correctly. Appropriate clinical training as well as a theoretical and practical proficiency of all the required operating techniques, including the use of this product, are prerequisites for the successful use of this product. The user is required to obtain information from the manufacturer if there is an unclear preoperative situation regarding the use of the product.

2.2.2 Sterility

The product is delivered in an unsterile condition.

- Clean the new product after removing its transport packaging and prior to its initial sterilization.

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*Note*  
Excessive mechanical stress may cause deformation and/or fracture.

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3.1 General safety information

*Note*  
Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

*Note*  
For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

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3.2 General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used. Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel. Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying. Additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at [www.a-k-i.org](http://www.a-k-i.org), link to "AKI-Brochures", "Red brochure".

3.3 Reusable products

Influences of the reprocessing which lead to damage to the product are not known. A careful visual and functional inspection before the next use is the best option to determine if a product is no longer functional, see Inspection.

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3.5 Cleaning/Disinfection

3.5.1 Product-specific safety information on the reprocessing method

Damage to or destruction of the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use cleaning agents and disinfectants approved for stainless steel in accordance with the manufacturer's instructions.
- Observe specifications regarding concentration, temperature and exposure time.
- Do not exceed the maximum allowable disinfection temperature of 95 °C.
- Use suitable cleaning/disinfecting agents if the product is disposed of in a wet condition. To prevent foaming and degradation of the efficacy of the process chemicals: prior to mechanical cleaning and disinfection, rinse the product thoroughly with running water.

3.5.2 Validated cleaning and disinfection procedure

| Validated procedure                                                                                 | Specific requirements                                                                                                                                                                                                                                                                                                      | Reference                                                                                                                                                                                                                                            |
|-----------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Manual cleaning with immersion disinfection                                                         | <ul style="list-style-type: none"><li>■ Suitable cleaning brush</li><li>■ 20 ml disposable syringe</li><li>■ Drying phase: Use lint free cloth or medical compressed air.</li></ul>                                                                                                                                        | Chapter Manual cleaning/disinfection and subsection: <ul style="list-style-type: none"><li>■ Chapter Manual cleaning with immersion disinfection</li></ul>                                                                                           |
| Manual pre-cleaning with brush and subsequent mechanical alkaline cleaning and thermal disinfection | <ul style="list-style-type: none"><li>■ Suitable cleaning brush</li><li>■ 20 ml disposable syringe</li><li>■ Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots).</li><li>■ Connect components with lumens and channels directly to the rinsing port of the injector carriage.</li></ul> | Chapter Mechanical cleaning/disinfection with manual pre-cleaning and subsection: <ul style="list-style-type: none"><li>■ Chapter Manual pre-cleaning with a brush</li><li>■ Chapter Mechanical alkaline cleaning and thermal disinfecting</li></ul> |

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3.6.1 Manual cleaning with immersion disinfection

| Phase | Step                  | T<br>[°C/°F] | t<br>[min] | Conc.<br>[%] | Water<br>quality | Chemical                                                       |
|-------|-----------------------|--------------|------------|--------------|------------------|----------------------------------------------------------------|
| I     | Disinfecting cleaning | RT<br>(cold) | >15        | 2            | D–W              | Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9* |
| II    | Intermediate rinse    | RT<br>(cold) | 1          | –            | D–W              | –                                                              |
| III   | Disinfection          | RT<br>(cold) | 5          | 2            | D–W              | Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9* |
| IV    | Final rinse           | RT<br>(cold) | 1          | –            | FD–W             | –                                                              |
| V     | Drying                | RT           | –          | –            | –                | –                                                              |

D–W: Drinking water  
FD–W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)  
RT: Room temperature  
\*Recommended: BBraun Stabimed fresh

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

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**Phase II**  
► Rinse/flush the product thoroughly (all accessible surfaces) under running water.  
► Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.  
► Drain any remaining water fully.

**Phase III**  
► Fully immerse the product in the disinfectant solution.  
► Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.  
► Rinse lumens at least 5 times at the beginning of the exposure time using an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

**Phase IV**  
► Rinse/flush the product thoroughly (all accessible surfaces).  
► Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.  
► Rinse lumens with an appropriate disposable syringe at least five times.  
► Drain any remaining water fully.

**Phase V**  
► Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

3.7 Mechanical cleaning/disinfection with manual pre-cleaning

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3.7.1 Manual pre-cleaning with a brush

| Phase | Step                  | T<br>[°C/°F] | t<br>[min] | Conc.<br>[%] | Water<br>quality | Chemical                                                       |
|-------|-----------------------|--------------|------------|--------------|------------------|----------------------------------------------------------------|
| I     | Disinfectant cleaning | RT<br>(cold) | >15        | 2            | D–W              | Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9* |
| II    | Rinsing               | RT<br>(cold) | 1          | –            | D–W              | –                                                              |

D–W: Drinking water  
RT: Room temperature  
\*Recommended: BBraun Stabimed fresh

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

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**Phase II**  
► Rinse/flush the product thoroughly (all accessible surfaces) under running water.  
► Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

3.7.2 Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

| Phase | Step                 | T<br>[°C/°F] | t<br>[min] | Water<br>quality | Chemical                                                                                                                                                                                                                                        |
|-------|----------------------|--------------|------------|------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| I     | Prerinse             | <25/77       | 3          | D–W              | –                                                                                                                                                                                                                                               |
| II    | Cleaning             | 55/131       | 10         | FD–W             | <div><div></div>Concentrate, alkaline:<ul style="list-style-type: none"><li>– pH ~ 13</li><li>– &lt;5 % anionic surfactant</li></ul><div><div></div>working solution 0.5%<ul style="list-style-type: none"><li>– pH = 11*</li></ul></div></div> |
| III   | Intermediate rinse   | >10/50       | 1          | FD–W             | –                                                                                                                                                                                                                                               |
| IV    | Thermal disinfecting | 90/194       | 5          | FD–W             | –                                                                                                                                                                                                                                               |
| V     | Drying               | –            | –          | –                | In accordance with the program for the cleaning and disinfecting machine                                                                                                                                                                        |

D–W: Drinking water  
FD–W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)  
\*Recommended: BBraun Helimatic Cleaner alkaline

- Check visible surfaces for residues after mechanical cleaning/disinfecting.

3.8 Inspection

- Allow the product to cool down to room temperature.
- Dry the product if it is wet or damp.

3.8.1 Visual inspection

- Ensure that all soiling has been removed. In particular, pay attention to mating surfaces, hinges, shafts, recessed areas, drill grooves and the sides of the teeth on rasps.
- If the product is dirty: repeat the cleaning and disinfection process.
- Check the product for damage, e.g. insulation or corroded, loose, bent, broken, cracked, worn or severely scratched and fractured components.
- Check the product for missing or faded labels.
- Check products having long, narrow geometries (in particular rotating instruments) for deformations.
- Check the product for damage to the spiral element.
- Check the cutting edges for continuity, sharpness, nicks and other damage.
- Check the surfaces for rough spots.
- Check the product for burrs that could damage tissue or surgical gloves.
- Check the product for loose or missing parts.
- Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical service.

3.8.2 Functional test

- Check that the product functions correctly.
- Check that all moving parts are working property (e.g. hinges, locks/latches, sliding parts etc.).
- Check for compatibility with associated products.
- Immediately put aside inoperative products and send them to Aesculap Technical Service, see Technical service.

3.9 Packaging

- Appropriately protect products with fine working tips.
- Place the product in its holder or on a suitable tray. Ensure that sharp edges are covered.
- Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- Ensure that the packaging provides sufficient protection against contamination of the product during storage.

3.10 Steam sterilization

- Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g., by opening any valves and faucets).
- Validated sterilization process
  - Steam sterilization using fractional vacuum process
  - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
  - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min
- If several devices are sterilized at the same time in the same steam sterilizer: Ensure that the maximum permitted load according to the manufacturers' specifications is not exceeded.

3.11 Storage

- Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

4. Technical service

**⚠ CAUTION**  
Modifications to medical equipment may result in the voiding of any guarantee/warranty claims and any approvals.  
► Do not modify the product.  
► For service and repairs, please contact your national B. Braun/Aesculap agency.

**Service address**  
Aesculap Technischer Service  
Am Aesculap-Platz  
78532 Tuttlingen / Germany  
Phone: +49 7461 95-1601  
Fax: +49 7461 16-2887  
E-Mail: ats@aesculap.de  
Other service addresses can be obtained from the address indicated above.

5. Disposal

**⚠ WARNING**  
Risk of infection due to contaminated products!  
► Adhere to national regulations when disposing of or recycling the product, its components and its packaging.  
  
**⚠ WARNING**  
Risk of injury from sharp and/or pointed products!  
► When disposing of or recycling the product, make sure the packaging prevents injuries from the product.

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|     |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
|-----|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| en  | <b>Instructions for use</b><br>Instruments with box-locks                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| USA | Note for U.S. users<br>This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at <a href="http://www.aesculapusaifus.com">www.aesculapusaifus.com</a> . If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost. |
| de  | <b>Gebrauchsanweisung</b><br>Instrumente mit Durchsteckschluss                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| fr  | <b>Mode d'emploi</b><br>Instruments à mécanisme de fermeture de sécurité                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| es  | <b>Instrucciones de manejo</b><br>Instrumental con tapa-pasador                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| it  | <b>Istruzioni per l'uso</b><br>Strumenti con meccanismo di incastro                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| pt  | <b>Instruções de utilização</b><br>Instrumentos com fecho de inserção                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| nl  | <b>Gebruiksaanwijzing</b><br>Instrumenten met insteekopening                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| da  | <b>Brugsanvisning</b><br>Gennemføringsværktøjer                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| nb  | <b>Bruksanvisning</b><br>Instrumenter med gjennomføringslås                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| sv  | <b>Bruksanvisning</b><br>Instrument med boxlås                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| fi  | <b>Käyttöohje</b><br>Läpivientien tiivisteillä varustetut instrumentit                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| et  | <b>Kasutusjuhend</b><br>Läbistava sulguriga instrumendid                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| lv  | <b>Lietošanas instrukcijas</b><br>Instrumenti ar punkcijas aizdari                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| lt  | <b>Naudojimo instrukcija</b><br>Instrumentai su perkisamuju užraktu                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| ru  | <b>Инструкция по применению</b><br>Инструменты со сквозным соединением замкового типа                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| cs  | <b>Návod k použití</b><br>Nástroje s prostrkávacím uzávěrem                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| pl  | <b>Instrukcja użytkowania</b><br>Instrumenty ze złączem przetykanym                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| sk  | <b>Návod na použitie</b><br>Nástroje s prestrkávacím uzáverom                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| hu  | <b>Használati útmutató</b><br>Műszerek átmenő csatlakozóval                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| sl  | <b>Navodila za uporabo</b><br>Instrumenti s priključkom za zapiranje                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| hr  | <b>Upute za uporabu</b><br>Instrumenti sa zatvaračem                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| ro  | <b>Manual de utilizare</b><br>Instrumente cu închidere prin întrepătrundere                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| bg  | <b>Употреба за употреба</b><br>Инструменти със затваряне чрез вмъкване                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| tr  | <b>Kullanım Kılavuzu</b><br>Geçme kilitli aletler                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| el  | <b>Οδηγίες χρήσης</b><br>Όργανα με ραβδωτή άκρη                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany  
Phone +49 (0) 7461 95-0 | Fax +49 (0) 7461 95-26 00 | [www.bbraun.com](http://www.bbraun.com)

AESCULAP® – a B. Braun brand

TA013729 2022-06 Change No. AE0061781

AESCULAP®  
Instruments with box-locks

1. About this document

Note

General risk factors associated with surgical procedures are not described in these instructions for use.

1.1 Scope

Multi-part instruments with one or more box-locks with and without industry-standard springs, which, except for the box-lock, have accessible and visible surfaces.

Note

The applicable CE mark for the product can be seen on the label or packaging of the product.

- For article specific instructions for use and material compatibility and lifetime information, see B. Braun eFU at eifu.bbraun.com

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⚠ WARNING

Indicates a possible threat of danger. If not avoided, minor or moderate injury may result.

⚠ CAUTION

Indicates a possible threat of material damage. If not avoided, the product may be damaged.

2. Clinical use

2.1 Areas of use and limitations of use

2.1.1 Intended use

The surgical instruments are intended for the universal use in various surgical disciplines.

2.1.2 Indications

Note

The manufacturer is not responsible for any use of the product against the specified indications and/or the described applications.

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Note

The user is obligated to report all severe events in connection with the product to the manufacturer and the responsible authorities of the state in which the user is located.

Notes on surgical procedures

It is the user's responsibility to ensure that the surgical procedure is performed correctly.

Appropriate clinical training as well as a theoretical and practical proficiency of all the required operating techniques, including the use of this product, are prerequisites for the successful use of this product.

The user is required to obtain information from the manufacturer if there is an unclear preoperative situation regarding the use of the product.

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⚠ WARNING

Risk of injury and/or malfunction!

- Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
- Always carry out a function test prior to each use of the product.

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3.1 General safety information

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note

Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.

Note

If there is no final sterilization, then a virucidal disinfectant must be used.

Note

For up-to-date information about reprocessing and material compatibility, see B. Braun eFU at eifu.bbraun.com

The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

3.2 General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying.

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3.3 Reusable products

Influences from reprocessing that lead to damage to the product are not known.

A careful visual and functional inspection before the next use provides the best opportunity to recognize a product that is no longer functional, see Inspection.

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- Use cleaning agents and disinfectants according to the manufacturer's instructions.
- Observe specifications regarding concentration, temperature and exposure time.
- Do not exceed the maximum allowable disinfection temperature of 95 °C.

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| Validated procedure                                   | Specific requirements                                                                                                                                                                                                                                                                                                                      | Reference                                                                                                                                                                 |
|-------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Manual cleaning with immersion disinfection           | <ul style="list-style-type: none"><li>■ Suitable cleaning brush</li><li>■ Disposable syringe 20 ml</li><li>■ When cleaning instruments with movable hinges, ensure that these are in an open position and, if applicable, move the hinge while cleaning.</li><li>■ Drying phase: Use a lint-free cloth or medical compressed air</li></ul> | Chapter Manual cleaning/disinfection and sub-section: <ul style="list-style-type: none"><li>■ Chapter Manual cleaning with immersion disinfection</li></ul>               |
| Mechanical alkaline cleaning and thermal disinfection | <ul style="list-style-type: none"><li>■ Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots).</li><li>■ Place the product on the tray with all product links and joints open.</li></ul>                                                                                                                   | Chapter Mechanical cleaning/disinfection and sub-section: <ul style="list-style-type: none"><li>■ Chapter Mechanical alkaline cleaning and thermal disinfection</li></ul> |

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| Phase | Step                  | T<br>[°C/°F] | t<br>[min] | Conc.<br>[%] | Water<br>quality | Chemical                                                       |
|-------|-----------------------|--------------|------------|--------------|------------------|----------------------------------------------------------------|
| I     | Disinfecting cleaning | RT<br>(cold) | >15        | 2            | D-W              | Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9* |
| II    | Intermediate rinse    | RT<br>(cold) | 1          | -            | D-W              | -                                                              |
| III   | Disinfection          | RT<br>(cold) | 5          | 2            | D-W              | Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9* |
| IV    | Final rinse           | RT<br>(cold) | 1          | -            | FD-W             | -                                                              |
| V     | Drying                | RT           | -          | -            | -                | -                                                              |

D-W:

Drinking water

FD-W:

Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT:

Room temperature

\*Recommended: BBraun Stabimed fresh

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3.7.1 Mechanical alkaline cleaning and thermal disinfection  
Machine type: single-chamber cleaning/disinfection device without ultrasound

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- \*Recommended: BBraun Helimatic Cleaner alkaline
- Check visible surfaces for residues after mechanical cleaning/disinfecting.

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3.8.1 Visual inspection

- Ensure that all soiling has been removed. In particular, pay attention to mating surfaces, hinges, shafts, recessed areas, drill grooves and the sides of the teeth on rasps.
- If the product is dirty: repeat the cleaning and disinfection process.
- Check the product for damage, e.g. insulation or corroded, loose, bent, broken, cracked, worn or severely scratched and fractured components.
- Check the product for missing or faded labels.
- Check the cutting edges for continuity, sharpness, nicks and other damage.
- Check the surfaces for rough spots.
- Check the product for burrs that could damage tissue or surgical gloves.
- Check the product for loose or missing parts.
- Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical service.

3.8.2 Functional test

- ⚠ CAUTION**  
Damage (metal cold welding/friction corrosion) to the product caused by insufficient lubrication!
- Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: STERILIT® I oil spray JG600 or STERILIT® I drip lubricator JG598).
  - Check that the product functions correctly.
  - Check that all moving parts are working property (e.g. hinges, locks/latches, sliding parts etc.).
  - Check for compatibility with associated products.
  - Immediately put aside inoperative products and send them to Aesculap Technical Service, see Technical service.

3.9 Packaging

- Appropriately protect products with fine working tips.
- Store products with ratchet locks fully opened or locked no further than in the first notch.
- Place the product in its holder or on a suitable tray. Ensure that sharp edges are covered.
- Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- Ensure that the packaging provides sufficient protection against contamination of the product during storage.

3.10 Steam sterilization

- Note*  
To avoid breakage due to stress crack corrosion, sterilize the instruments with the lock fully open or locked no further than on the first ratchet tooth.
- Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g., by opening any valves and faucets).
  - Validated sterilization process
    - Steam sterilization in fractionated vacuum process
    - Steam sterilizer in accordance with DIN EN 285 and validated in accordance with DIN EN ISO 17665
    - Sterilization in fractionated vacuum process at 134 °C, holding time 5 min.
  - If several devices are sterilized at the same time in the same steam sterilizer: Ensure that the maximum permitted load according to the manufacturers' specifications is not exceeded.

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4. Technical service

- ⚠ CAUTION**  
Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.
- Do not modify the product.
  - For service and repairs, please contact your national B. Braun/Aesculap agency.

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5. Disposal

- ⚠ WARNING**  
Risk of infection due to contaminated products!
- Adhere to national regulations when disposing of or recycling the product, its components and its packaging.

- ⚠ WARNING**  
Risk of injury due to sharp-edged and/or pointed products!
- When disposing of or recycling the product, ensure that the packaging prevents injury by the product.

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TA013729    2022-06    Change No. AE0061781

|     |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
|-----|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| en  | <b>Instructions for use/Technical description</b><br>Ring-type scissors and spring-type scissors                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| USA | Note for U.S. users<br>This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at <a href="http://www.aesculapusaifus.com">www.aesculapusaifus.com</a> . If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost. |
| de  | <b>Gebrauchsanweisung/Technische Beschreibung</b><br>Ring- und Federscheren                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| fr  | <b>Mode d'emploi/Description technique</b><br>Ciseaux à anneau et à ressort                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| es  | <b>Instrucciones de manejo/Descripción técnica</b><br>Tijeras anulares y elásticas                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| it  | <b>Istruzioni per l'uso/Descrizione tecnica</b><br>Forbici ad anello e a molla                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| pt  | <b>Instruções de utilização/Descrição técnica</b><br>Tesouras com anéis e de mola                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| nl  | <b>Gebruiksaanwijzing/Technische beschrijving</b><br>Ring- en veerscharen                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| da  | <b>Brugsanvisning/Teknisk beskrivelse</b><br>Ring- og fjedersakse                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| sv  | <b>Bruksanvisning/Teknisk beskrivning</b><br>Ring- och fjädersaxar                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| fi  | <b>Käyttöohje/Tekninen kuvaus</b><br>Rengas- ja jousisakset                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| lv  | <b>Lietošanas instrukcijas/tehniskais apraksts</b><br>Šķēres ar apaļiem rokturiem un atsperšķēres                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| lt  | <b>Naudojimo instrukcija/techninis aprašas</b><br>Žiedinės ir spyruoklinės žirkės                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| ru  | <b>Инструкция по применению/Техническое описание</b><br>Изогнутые и пружинные ножницы                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| cs  | <b>Návod k použití/Technický popis</b><br>Kroužkové a pružinové nůžky                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| pl  | <b>Instrukcja użytkowania/Opis techniczny</b><br>Nożyczki chirurgiczne i mikrochirurgiczne                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| sk  | <b>Návod na použitie/Technický opis</b><br>Nožnice s okami a pružinové nožnice                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| hu  | <b>Használati útmutató/Műszaki leírás</b><br>Gyűrűs és rugós ollók                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| sl  | <b>Navodila za uporabo/Tehnični opis</b><br>Škarje z okroglim ročajem in vzmetne škarje                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| hr  | <b>Upute za uporabu/Tehnički opis</b><br>Škare s prstenom i oprugom                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| ro  | <b>Manual de utilizare/Descriere tehnică</b><br>Foarfeci cu inel și arc                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| bg  | <b>Упътване за употреба/Техническо описание</b><br>Пръстеновидни и пружинни ножници                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| tr  | <b>Kullanım Kılavuzu/Teknik açıklama</b><br>Halka ve yaylı makaslar                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| el  | <b>Οδηγίες χρήσης/Τεχνική περιγραφή</b><br>Ψαλίδι με δακτυλίδι και ψαλίδι ελατηρίου                                                                                                                                                                                                                                                                                                                                                                                                                                                          |

## B | BRAUN

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany  
Phone +49 (0) 7461 95-0 | Fax +49 (0) 7461 95-26 00 | [www.bbraun.com](http://www.bbraun.com)

AESCULAP® – a B. Braun brand

TA015820 2021-01 Change No. 63662

Ring-type scissors and spring-type scissors

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1.1 Scope

These instructions for use apply to ring-type and spring-type scissors in nearly all surgical disciplines.

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2.1 Areas of use and limitations of use

2.1.1 Intended use

**Surgical scissors**  
The instruments are used to cut tissue and/or medical materials and supplies.

**Dissecting scissors**  
The instruments are used to cut and/or dissect tissue.

**Nail scissors**  
The instruments are used to cut or split finger nails and toe nails and/or cuticles.

**Bandage scissors and material scissors**  
The instruments are used to cut medical materials and supplies and/or clothing.

**Micro scissors**  
The instruments are used to cut and/or dissect tissue during micro surgical procedures.

**2.1.2 Indications**  
The instruments are used in a multitude of surgical procedures and interventions in almost all surgical disciplines, see Intended use.

**2.1.3 Contraindications**  
No contraindications for the product are currently known.

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**2.2.2 Sterility**  
The product is delivered in an unsterile condition.  
► Clean the new product after removing its transport packaging and prior to its initial sterilization.

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3.5 Cleaning/Disinfection

**3.5.1 Product-specific safety information on the reprocessing method**  
Damage to or destruction of the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!  
► Use cleaning agents and disinfectants according to the manufacturer's instructions.  
► Observe specifications regarding concentration, temperature and exposure time.  
► Do not exceed a disinfection temperature of 96 °C.

For products with plasma layers (e. g. Noir instruments), the layer can be damaged or worn if special cleaning procedures with oxidising chemicals (e. g. hydrogen peroxide H<sub>2</sub>O<sub>2</sub>) are used.  
► Do not use oxidizing chemicals for cleaning.  
► For wet disposal, use suitable cleaning agents/disinfectants. To prevent foam formation and reduced effectiveness of the process chemicals: Prior to mechanical cleaning and disinfection, rinse the product thoroughly under running water

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| Validated procedure                                   | Specific requirements                                                                                                                                                                                                                                                                                | Reference                                                                                                                   |
|-------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------|
| Manual cleaning with immersion disinfection           | <ul style="list-style-type: none"><li>► Use a suitable cleaning brush.</li><li>► Keep working ends open when cleaning.</li><li>► Clean products with movable hinges in the open position or while moving the joints.</li><li>► Drying phase: Use lint free cloth or medical compressed air</li></ul> | Chapter Manual cleaning/disinfection and subsection:<br>■ Chapter Manual cleaning with immersion disinfection               |
| Mechanical alkaline cleaning and thermal disinfection | <ul style="list-style-type: none"><li>► Place product on screen basket suitable for cleaning (make sure all areas will be reached by water jets).</li><li>► Keep working ends open when cleaning.</li><li>► Placing the product on the sterilization tray with the hinge open.</li></ul>             | Chapter Mechanical cleaning/disinfection and subsection:<br>■ Chapter Mechanical alkaline cleaning and thermal disinfection |

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| Phase | Step                  | T<br>[°C/°F] | t<br>[min] | Conc.<br>[%] | Water<br>quality | Chemical                                                       |
|-------|-----------------------|--------------|------------|--------------|------------------|----------------------------------------------------------------|
| I     | Disinfecting cleaning | RT<br>(cold) | >15        | 2            | D–W              | Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9* |
| II    | Intermediate rinse    | RT<br>(cold) | 1          | –            | D–W              | –                                                              |
| III   | Disinfection          | RT<br>(cold) | 5          | 2            | D–W              | Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9* |
| IV    | Final rinse           | RT<br>(cold) | 1          | –            | FD-W             | –                                                              |
| V     | Drying                | RT           | –          | –            | –                | –                                                              |

D–W: Drinking water  
FD–W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)  
RT: Room temperature  
\*Recommended: BBraun Stabimed fresh  
► Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

- Phase I**
- Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
  - Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
  - If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
  - Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
  - Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

- Phase II**
- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
  - Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
  - Drain any remaining water fully.

- Phase III**
- Fully immerse the product in the disinfectant solution.
  - Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
  - Rinse lumens at least 5 times at the beginning of the exposure time using an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

- Phase IV**
- Rinse/flush the product thoroughly (all accessible surfaces).
  - Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
  - Rinse lumens with an appropriate disposable syringe at least five times.
  - Drain any remaining water fully.

- Phase V**
- Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

3.7 Mechanical cleaning/disinfection

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3.7.1 Mechanical alkaline cleaning and thermal disinfection

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| Phase | Step                 | T<br>[°C/°F] | t<br>[min] | Water<br>quality | Chemical/Note                                                                                                                                                                                                                                    |
|-------|----------------------|--------------|------------|------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| I     | Pre-rinse            | <25/77       | 3          | D-W              | -                                                                                                                                                                                                                                                |
| II    | Cleaning             | 55/131       | 10         | FD-W             | <div><div></div>Concentrate, alkaline:<ul style="list-style-type: none"><li>- pH = 13</li><li>- &lt;5 % anionic surfactant</li></ul><div><div></div>0.5 % working solution<ul style="list-style-type: none"><li>- pH = 11*</li></ul></div></div> |
| III   | Intermediate rinse   | >10/50       | 1          | FD-W             | -                                                                                                                                                                                                                                                |
| IV    | Thermal disinfecting | 90/194       | 5          | FD-W             | -                                                                                                                                                                                                                                                |
| V     | Drying               | -            | -          | -                | According to the program for cleaning and disinfection device                                                                                                                                                                                    |

D-W: Drinking water  
FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)  
\*Recommended: BBraun Helimatic Cleaner alkaline

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3.8.1 Visual inspection

- Make sure all dirt has been removed. In particular, pay attention to, for example, mating surfaces, hinges, shafts, recessed areas and drilled grooves.
- If the product is dirty: repeat the cleaning and disinfection process.
- Check product for damage, e.g., corroded, loose, bent, broken, cracked, worn, heavily scratched or missing parts.
- Check the product for missing or faded labels.
- Check the cutting edges for continuity, sharpness, nicks and other damage.
- Check the surfaces for rough spots.
- Check the product for burrs that could damage tissue or surgical gloves.
- Check the product for loose or missing parts.
- Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical service.

3.8.2 Functional test

 CAUTION

Damage (metal cold welding/friction corrosion) to the product caused by insufficient lubrication!

- Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: STERILIT® I oil spray JG600 or STERILIT® I drip lubricator JG598).
- Check that the product functions correctly.
- Check that all moving parts are working properly (e.g. hinges, locks/latches, sliding parts etc.).
- Make sure that the scissor blades provide light resistance when closing.
- Immediately put aside inoperative products and send them to Aesculap Technical Service, see Technical service.

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3.10 Steam sterilization

- Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g., by opening any valves and faucets).
- Validated sterilization process
  - Steam sterilization in fractionated vacuum process
  - Steam sterilizer in accordance with DIN EN 285 and validated in accordance with DIN EN ISO 17665
  - Sterilization in fractionated vacuum process at 134 °C, holding time 5 min
- If several devices are sterilized at the same time in the same steam sterilizer: Ensure that the maximum permitted load according to the manufacturers' specifications is not exceeded.

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4. Technical service

 CAUTION

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

- Do not modify the product.
- For service and repairs, please contact your national B. Braun/Aesculap agency.

Service addresses

Aesculap Technischer Service

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78532 Tuttlingen / Germany

Phone: +49 7461 95-1601

Fax: +49 7461 16-2887

E-Mail: ats@aesculap.de

Other service addresses can be obtained from the address indicated above.

5. Disposal

 WARNING

Risk of infection due to contaminated products!

- Adhere to national regulations when disposing of or recycling the product, its components and its packaging.

 WARNING

Risk of injury due to sharp-edged and/or pointed products!

- When disposing of or recycling the product, ensure that the packaging prevents injury by the product.

*Note*

The user institution is obliged to reprocess the product before its disposal, see Validated reprocessing procedure.

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Aesculap®

Single-piece instruments made from high-grade steel, tin and plastics suitable for alkaline reprocessing

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1.1 Scope

Single-piece instruments, of which all surfaces are directly accessible and visible. The instruments do not contain aluminum or nickel-silver components and are not nickel-plated or chrome-plated.

*Note*  
The applicable CE marking for the product can be found on the label or packaging of the product.

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2. Clinical use

2.1 Areas of use and limitations of use

**2.1.1 Intended use**  
The surgical instruments are intended for the universal use in various surgical disciplines.

**2.1.2 Indications**  
*Note*  
The manufacturer is not responsible for any use of the product against the specified indications and/or the described applications.

For indications, see Intended use.

**2.1.3 Absolute contraindications**  
No known absolute contraindications.

**2.1.4 Relative contraindications**  
The following conditions, individual or combined, can lead to delayed healing or compromise the success of the operation:  
■ Medical or surgical conditions (e.g. comorbidities) which could hinder the success of the operation.  
In the presence of relative contraindications, the user decides individually regarding the use of the product.

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**2.2.2 Sterility**  
The product is delivered in an unsterile condition.  
► Clean the new product after removal of the transport packaging (including the protective cover of the working tips) and before first sterilization.

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3.1 General safety information

*Note*  
Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.  
*Note*  
For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

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3.5 Cleaning/Disinfection

**3.5.1 Product-specific safety information on the reprocessing method**  
Damage to or destruction of the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!  
► Use cleaning and disinfecting agents approved for high-grade steel and plastics according to manufacturer's instructions.  
► Observe specifications regarding concentration, temperature and exposure time.  
► Do not exceed the maximum allowable disinfection temperature of 95 °C.  
► If there are bone, tissue or ancillary material (e.g. plaster, bone cement) residues: pre-clean the product manually (using a cleaning brush).

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| Validated procedure                                   | Specific requirements                                                                                                      | Reference                                                                                                                   |
|-------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------|
| Manual cleaning with immersion disinfection           | ■ Suitable cleaning brush<br>■ Disposable syringe 20 ml<br>■ Drying phase: Use a lint-free cloth or medical compressed air | Chapter Manual cleaning/disinfection and subsection:<br>■ Chapter Manual cleaning with immersion disinfection               |
| Mechanical alkaline cleaning and thermal disinfection | ■ Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots).                                   | Chapter Mechanical cleaning/disinfection and subsection:<br>■ Chapter Mechanical alkaline cleaning and thermal disinfection |

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| Phase         | Step                                                                                                        | T [°C/°F] | t [min] | Conc. [%] | Water quality | Chemical                                                       |
|---------------|-------------------------------------------------------------------------------------------------------------|-----------|---------|-----------|---------------|----------------------------------------------------------------|
| I             | Disinfecting cleaning                                                                                       | RT (cold) | >15     | 2         | D-W           | Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9* |
| II            | Intermediate rinse                                                                                          | RT (cold) | 1       | -         | D-W           | -                                                              |
| III           | Disinfection                                                                                                | RT (cold) | 5       | 2         | D-W           | Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9* |
| IV            | Final rinse                                                                                                 | RT (cold) | 1       | -         | FD-W          | -                                                              |
| V             | Drying                                                                                                      | RT        | -       | -         | -             | -                                                              |
| D-W:          | Drinking water                                                                                              |           |         |           |               |                                                                |
| FD-W:         | Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least) |           |         |           |               |                                                                |
| RT:           | Room temperature                                                                                            |           |         |           |               |                                                                |
| *Recommended: | BBraun Stabimed fresh                                                                                       |           |         |           |               |                                                                |

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3.10 Steam sterilization

- Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g., by opening any valves and faucets).
- Validated sterilization process
  - Steam sterilization in fractionated vacuum process
  - Steam sterilizer in accordance with DIN EN 285 and validated in accordance with DIN EN ISO 17665
  - Sterilization in fractionated vacuum process at 134°C, holding time 5 min
- If several devices are sterilized at the same time in the same steam sterilizer: Ensure that the maximum permitted load according to the manufacturers' specifications is not exceeded.

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4. Technical service

⚠ CAUTION

- Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.
- Do not modify the product.
  - For service and repairs, please contact your national B. Braun/Aesculap agency.

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5. Disposal

⚠ WARNING

- Risk of infection due to contaminated products!
- Adhere to national regulations when disposing of or recycling the product, its components and its packaging.

⚠ WARNING

- Risk of injury due to sharp-edged and/or pointed products!
- When disposing of or recycling the product, ensure that the packaging prevents injury by the product.

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